

 **FX-V**  
Fowler Bed



Vertical  
Siderails

Fowler Position

Smooth  
Crank System



*#polymerheadandfootpanels #standardconfiguration*



# Fowler Bed

## Clinical Positions



Backrest



Kneerest



DVT Platform

## Features

- Backrest and kneerest adjustable on crank mechanism.
- Four section perforated mattress platform.
- The adjustable calf rest adapts to obtain a good DVT platform.
- The Fowler bed allows the caregiver to obtain basic clinical positions.
- Anti-microbial powder coating.

## Specifications

Model	FX-V
External dimension (L X W) mm	2200 X 955
External dimension with Vertica siderails (L X W) mm	2200 X 1000
Mattress platform dimension (L X W) mm	1980 X 805
Fixed height (H) mm (without mattress)	600
Backrest angle	0° to 65°
Kneerest angle	0° to 40°
SWL - Safe working load	135 kgs
Patient weight bearing capacity	105 kgs

Engineering variation ±2°, ±10mm

## Standard Accessories

Polymer moulded head and foot panels (PP3) - easily removable
Four section 40 density foam mattress of 100 mm thickness covered with rexine (M3)
Removable chromium plated forged handle
IV pole provision
Mattress arrestor

## Optional Accessories

3/4 Length Vertica siderails
125 mm Ø high grade synthetic body casters with brake on two diagonally opposite casters (C3)
SS 304 telescopic IV pole with two hooks
Removable polymer moulded handle
Urine bag holder

Subject to technical changes without notice. Accessories in the picture may not be part of the standard configuration. The product colours in the pictures may differ slightly from the actual product.

### MUMBAI OFFICE

AMSON BIOMED Pvt. Ltd.,

1st Floor, Lal Bungalow, Upashray Lane, Ghatkopar East, Mumbai 400077

Phone +91 9820032799 | info@amsonbiomed.com | www.amsonbiomed.com

