

Iris-P

Motorized Examination Table



Easy Access

Height Adjustment
by Foot Switch Control

Low Height for
Patient Comfort



#dignifiedpatientexamination

Iris-P

Motorized Examination Table

Features

- Motorized height adjustment on foot switch control.
- The low height ensures easy access to the patient for sitting and alighting from the exam table.
- Backrest adjustable on gas spring.
- Storage drawer accessible from both sides.
- Adjustable tray for BP apparatus.
- Integrated paper roll holder.
- Anti-microbial powder coating.

Specifications

Model	Iris-P
External dimension (L X W) mm	1980 X 780
Mattress platform dimension (L X W) mm	1980 X 650
Exam table height adjustment with upholstery (H) mm	535 - 955
Headrest angle	0° to 60°
SWL - Safe working load	170 kgs
Patient weight bearing capacity	150 kgs

Engineering variation $\pm 2^\circ$, ± 10 mm

Electrical Parameters

Liquid ingress protection	IP54
Shock protection	Class I type B
Power in	Nominal 230V AC Max 2A 90V-300V

Safety Standards

IEC 60601-2-52



Subject to technical changes without notice. Accessories in the picture may not be part of the standard configuration. The product colours in the pictures may differ slightly from the actual product.