

➤ Patient Care Solutions



ST43A

ST43A SS Mobile IV Stand with two hooks

Stainless steel tubular construction mounted on five pronged rectangular tubular base
50mm (2") Ø twin wheel synthetic body casters
Stainless steel rod with two hooks
Height adjustment (H): 1315mm to 2440mm (51.5" to 96")

MT32 SS Mayo's Table

Stainless steel tray height adjustment: 890 mm to 1550 mm (35" to 61")
Tray dimension (LXW) : 560 mm X 400 mm (22" X 16")
Stainless steel tubular frame mounted on 50mm (2") Ø twin wheel synthetic body casters



MT32



TTB11

TTB11 Two Tier Bowl Stand

Overall approximate height (H): 850 mm (33")
Mild steel tubular frame mounted on three 50mm (2") Ø twin wheel synthetic body casters
Two stainless steel bowls of approximately 355mm (14") Ø

➤ Patient Care Solutions

SS21 Stainless Steel Stool

Overall approximate dimension (LXWXH):
305mm X 305mm X 505mm (12" X 12" X 20")

Base area dimension (LXW): 410 mm X 410 mm (16" X 16")

Stainless steel structure



SS21



RS02

RS02 Stainless Steel Revolving Stool

Height adjustment: 460 mm to 665 mm (18" to 26")

Stainless steel structure with revolving top

FD14 Stainless Steel Foldable food table

Overall approximate dimension (LXWXH):
560 mm X 400 mm X 260 mm (22" X 16" X 10.5")

Stainless steel tray



FD14



KB12

KB12 Stainless Steel Kick Bucket

Stainless steel bowl 355 mm (14") Ø

The tubular framework is fitted with three 125mm (5") Ø twin wheel synthetic body casters

Overall height: 335 mm (15")

MUMBAI OFFICE

AMSON BIOMED Pvt. Ltd.,

1st Floor, Lal Bungalow, Upashray Lane, Ghatkopar East, Mumbai 400077

Phone +91 9820032799 | info@amsonbiomed.com | www.amsonbiomed.com

