

Revero-G

Bedside Table



Features

- Revero with two section membrane top. The bigger section is inclinable on a pivot/ratchet as desired and the smaller section is fixed.
- Height adjustment by Gas Spring Mechanism gives the patient comfortable table levels for eating and resting. In addition, smaller medical equipment and any routine testing apparatus can be placed on the table top.
- Raised edges to stop objects from falling off the table top.
- Anti-microbial powder coating.

Specifications

Model	Revero-G
External dimension (L X W) mm	960 X 410
Height adjustment (H) mm	785 to 1065
Inclinable top section dimension (L X W) mm	615 X 355
Fixed top section dimension (L X W) mm	205 X 355
Tilting angle	0° to 45°
Casters	50 mm Ø twin wheel synthetic body casters

Engineering variation ±2°, ±10mm

#easyheightadjustment

Revero-S

Bedside Table



Features

- Revero with two section membrane top. The bigger section is inclinable on a pivot/ratchet as desired and the smaller section is fixed.
- Height adjustment by Smooth Screw Mechanism gives the patient comfortable table levels for eating and resting. In addition, smaller medical equipment and any routine testing apparatus can be kept.
- Raised edges to stop objects from falling off the table top.
- Anti-microbial powder coating.

Specifications

Model	Revero-S
External dimension (L X W) mm	960 X 410
Height adjustment (H) mm	715 to 1080
Inclinable top section dimension (L X W) mm	615 X 355
Fixed top section dimension (L X W) mm	205 X 355
Tilting angle	0° to 45°
Casters	50 mm Ø twin wheel synthetic body casters

Engineering variation $\pm 2^\circ$, $\pm 10\text{mm}$

#easyheightadjustment