

 **SFX-M**  
Semi-Fowler Bed



  
3/4 Length  
Siderails

  
Raised  
Backrest

  
Smooth  
Crank System



*#polymerheadandfootpanels #standardconfiguration*

# SFX-M

## Semi-Fowler Bed

### Clinical Positions



### Features

- Backrest adjustable on crank mechanism.
- Two section perforated mattress platform.
- Anti-microbial powder coating.

### Specifications

Model	SFX-M
External dimension (L X W) mm	2200 X 955
External dimension with 3/4 length siderails (L X W) mm	2200 X 1000
Mattress platform dimension (L X W) mm	1980 X 805
Fixed height (H) mm (without mattress)	600
Backrest angle	0° to 65°
SWL - Safe working load	135 kgs
Patient weight bearing capacity	105 kgs

Engineering variation  $\pm 2^\circ$ ,  $\pm 10$ mm

### Standard Accessories

Polymer moulded head and foot panels (PP3) - easily removable
Two section 40 density foam mattress of 100 mm thickness covered with rexine (M4)
Removable chromium plated forged handle
IV pole provision
Mattress arrestor

### Optional Accessories

3/4 Length siderails
125 mm $\varnothing$ high grade synthetic body casters with brake on two diagonally opposite casters (C3)
SS 304 telescopic IV pole with two hooks
Removable polymer moulded handle
Urine bag holder

Subject to technical changes without notice. Accessories in the picture may not be part of the standard configuration. The product colours in the pictures may differ slightly from the actual product.