


# SS Crash Cart-S



**The Stainless Steel Crash Cart is designed to carry life saving equipment and drugs**

*#emergencycare*

# SS Crash Cart-S

## Specifications

Model	SS Crash Cart-S
External dimension (L X W) mm	950 X 465
Platform dimension (L X W) mm	630 X 410
Height (H) mm	1650
Defibrillator tray dimension (L X W) mm	650 X 195
Cardiac board dimension (L X W) mm	670 X 350
Shelf dimension (L X W) mm	670 X 410

Engineering variation ±10mm

## Standard Accessories

Three shelves
Two storage plastic boxes with three drawers each
12 Storage plastic bins
CPR board cum extension table
SS 304 IV pole
SS 304 protective siderails
Oxygen tank holder
125 mm Ø synthetic body casters with brake on two diagonally opposite casters (C5)

Subject to technical changes without notice. Accessories in the picture may not be part of the standard configuration. The product colours in the pictures may differ slightly from the actual product.

### MUMBAI OFFICE

AMSON BIOMED Pvt. Ltd.,

1st Floor, Lal Bungalow, Upashray Lane, Ghatkopar East, Mumbai 400077

Phone +91 9820032799 | info@amsonbiomed.com | www.amsonbiomed.com

