

SS Trolleys

Dressing & Instrument Trolleys



IT24



IT68



DT67

Specifications	Stainless Steel Instrument Trolley		
	IT24	IT68	DT67
Model	IT24	IT68	DT67
External approximate dimension (L X W X H) inches	27 X 18 X 36	30 X 20 X 36	40 X 20 X 36
External dimension (LXWXH) mm	680 X 450 X 920	760 X 510 X 920	1020 X 510 X 920
Shelf dimension (LXW)	680 X 450	755 X 505	755 X 505
Shelves	2 nos	2 nos	2 nos
Safety siderails	3 sides	4 sides	4 sides
Casters	125 mm Ø synthetic body casters with brake on two diagonally opposite casters (C5)		

Engineering variation ±10mm

Subject to technical changes without notice. Accessories in the picture may not be part of the standard configuration. The product colours in the pictures may differ slightly from the actual product.

SS Trolleys

Instrument Trolleys



IT25



IT25A



IT25B

Specifications	Stainless Steel Instrument Trolley		
	IT25	IT25A	IT25B
Model	IT25	IT25A	IT25B
External approximate dimension (L X W X H) inches	27 X 18 X 36	42 X 24 X 36	36 X 20 X 36
External dimension (LXWXH) mm	680 X 450 X 920	1070 X 610 X 920	915 X 510 X 920
Shelf dimension (LXW)	680 X 450	1070 X 610	915 X 510
Shelves	2 nos	2 nos	2 nos
Safety siderails	3 sides	3 sides	3 sides
Casters	125 mm Ø synthetic body casters with brake on two diagonally opposite casters (C5)		

Engineering variation ±10mm

Subject to technical changes without notice. Accessories in the picture may not be part of the standard configuration. The product colours in the pictures may differ slightly from the actual product.