

Spacio-M

Bedside Cabinet



Features

- Cabinet top with rounded corners. The front side of the drawer and door in membrane finish.
- Three sided SS rail on the top to prevent objects from falling down.
- SS handles for the drawer and the door.
- Protective bumpers at back of the cabinet.
- Anti-microbial powder coating.

Specifications

Model	Spacio-M
External dimension (L X W X H) mm	450 X 400 X 870
External dimension with document holder (L X W X H) mm	495 X 465 X 870
Drawer dimension (H X W X D) mm	90 X 335 X 370
Cabinet dimension (H X W X D) mm	375 X 385 X 400
Casters	50 mm Ø twin wheel synthetic body casters

Engineering variation ±10mm

Optional accessories

Document holder

#comfortablestoragespace

Spacio-S

Bedside Cabinet



Features

- Cabinet top made of SS 304, drawer and door made of MS.
- Three sided raised edges on the top to prevent objects from falling down.
- Protective bumpers at back of the cabinet.
- Anti-microbial powder coating.

Specifications

Model	Spacio-S
External dimension (L X W X H) mm	400 X 425 X 845
Drawer dimension (H X W X D) mm	100 X 330 X 390
Cabinet dimension (H X W X D) mm	370 X 360 X 400
Casters	50 mm Ø twin wheel synthetic body casters

Engineering variation ±10mm

#compactstoragespace