

➤ Stretcher Trolley-ST1



Features

- Removable two section stretcher top made of HPL (high pressure laminate).
- Backrest adjustment on ratchet.
- Ergonomically designed handles for better manoeuvring.
- Anti-microbial powder coating.

Specifications

Model	ST1
External dimension (L X W) mm	1920 X 650
Stretcher dimension (L X W) mm	1785 X 555
Fixed stretcher height (H) mm (without mattress)	840
Backrest angle	65°
SWL - Safe working load	135 kgs
Patient weight bearing capacity	105 kgs

Engineering variation $\pm 2^\circ$, ± 10 mm

Standard Accessories

200 mm \varnothing casters with brake on two diagonally opposite casters (C6)
Mattress of 50 mm foam with rexine
Safety belts
MS swing away type siderails
SS 304 IV pole - fixed height
MS oxygen cylinder cage

Optional Accessories

MS storage tray

Stretcher Trolley-ST2



Features

- Removable aluminium stretcher top.
- Ergonomically designed handles for better manoeuvring.
- Anti-microbial powder coating.

Specifications

Model	ST2
External dimension (L X W) mm	2070 X 650
Stretcher dimension (L X W) mm	1820 X 560
Fixed stretcher height (H) mm (without mattress)	830
SWL - Safe working load	135 kgs
Patient weight bearing capacity	130 kgs

Engineering variation ± 10 mm

Standard Accessories

200 mm \varnothing casters with brake on two diagonally opposite casters (C6)
Safety belts

Optional Accessories

Mattress of 50 mm foam with rexine
MS swing away type siderails
SS 304 IV pole - fixed height
MS oxygen cylinder cage
MS storage tray

#safeexpress

MUMBAI OFFICE

AMSON BIOMED Pvt. Ltd.,

1st Floor, Lal Bungalow, Upashray Lane, Ghatkopar East, Mumbai 400077

Phone +91 9820032799 | info@amsonbiomed.com | www.amsonbiomed.com

